

Production Order Form – Mould Detect for Encanta

I. Customer details

Order Order no _____

Rush Order

Date Ordered date _____
Requested delivery date _____

End user First name _____
Surname _____

Hospital / Clinic Phone _____
Hospital / Clinic _____
Address _____
City _____ Country _____ Post code _____
E-mail _____


End user info Phone _____ Birth date _____
Address _____
City _____ Country _____ Post code _____
E-mail _____

Audiogram

	125 Hz	250 Hz	500 Hz	1000 Hz	2000 Hz	4000 Hz	8000 Hz
Left							
Right							

2. Product configuration

Encanta miniRITE Micro Mould



	Retention Style																					
	Micro		Canal		Double Canal Lock		Canal Lock		Half Shell		Full Shell		Half Skeleton		Skeleton		Skeleton Wing Lock		Helix Lock		Double Canal Lock Short	
Earpiece material	L	R	L	R	L	R	L	R	L	R	L	R	L	R	L	R	L	R	L	R	L	R
Hard acrylic											N/A											
OtoTherm, shore 70																						
Soft silicone, shore 40					N/A						N/A						N/A	N/A			N/A	
Soft silicone, shore 60					N/A						N/A						N/A	N/A			N/A	

Mandatory Options

Ventilation
 Left: 0.8 (S) 1.4 (M) 2.4 (L) 3.0 (XL) No vent
 Right: 0.8 (S) 1.4 (M) 2.4 (L) 3.0 (XL) No vent

Fitting Level
 Left: 60 85
 Right: 60 85

Micro Mould Type
 Left: Short Long
 Right: Short Long


















Optional

Speaker Unit Length
 Left: 0 1 2 3 4 5
 Right: 0 1 2 3 4 5

Production Order Form – Mould Detect

Encanta miniRITE

LiteTip

	Retention Style																
	Micro		Canal		Canal Lock		Half Skeleton		Skeleton		Double Canal Lock		Skeleton Wing Lock		Double Canal Lock Short		
																	
Earpiece material		L	R	L	R	L	R	L	R	L	R	L	R	L	R	L	R
Hard acrylic																	
OtoTherm, shore 70																	

Mandatory Options




















Ventilation
 Left: 0.7 (S) 0.8 (M) 1.2 (L) 1.5 (XL) No vent
 Right: 0.7 (S) 0.8 (M) 1.2 (L) 1.5 (XL) No vent

Fitting Level
 Left: 60 85
 Right: 60 85

Optional

Speaker Unit Length
 Left: 0 1 2 3 4 5
 Right: 0 1 2 3 4 5

MicroShell

	Retention Style																					
	Micro ⁽¹⁾		Canal		Double Canal Lock		Canal Lock		Half Skeleton		Skeleton		Skeleton Wing Lock		Helix Lock		Double Canal Lock Short		Half Shell ⁽²⁾		Full Shell ⁽²⁾	
																						
Earpiece material		L	R	L	R	L	R	L	R	L	R	L	R	L	R	L	R	L	R	L	R	
Hard acrylic																						

Mandatory Options

Ventilation Cast (Round)
 Left: 0.8 1.4 1.8 2.4 3.0⁽³⁾ Open vent⁽³⁾ No vent
 Right: 0.8 1.4 1.8 2.4 3.0⁽³⁾ Open vent⁽³⁾ No vent

Ventilation Collection - Only MiniFit Detect 100 & 105 FL
 Left: 1.0 1.4 2.4 3.0⁽³⁾
 Right: 1.0 1.4 2.4 3.0⁽³⁾

MicroShell Type
 Left: Short⁽⁴⁾ Long
 Right: Short⁽⁴⁾ Long

Fitting Level
 Left: 60 85 100 105
 Right: 60 85 100 105

Fixed Wire Length
 Left: 1 2 3 4 5
 Right: 1 2 3 4 5

¹ Micro retention style not available for FL 100 & FL 105
² Half Shell & Full Shell available only for FL100 & FL105
³ Not available for FL105
⁴ Not available for FL100 & FL105

3. Additional information

Laser Marking

Standard (Serial number)
 Standard & Initials: _____
 Other: _____

Colour

		Colour																															
		Transparent		Beige ⁽¹⁾		Light Brown ⁽²⁾		Medium Brown ⁽²⁾		Dark Brown ⁽²⁾		Transparent Blue		Transparent Red		Rose		Blue		Red		Orange		Yellow		Green		Pink		White		Black ⁽³⁾	
Earpiece material		L	R	L	R	L	R	L	R	L	R	L	R	L	R	L	R	L	R	L	R	L	R	L	R	L	R	L	R	L	R		
Hard acrylic		<input type="checkbox"/>	<input type="checkbox"/>															N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A				
Soft silicone, shore 40		<input type="checkbox"/>	<input type="checkbox"/>			N/A	N/A	N/A	N/A	N/A	N/A																	N/A					
Soft silicone, shore 60		<input type="checkbox"/>	<input type="checkbox"/>	N/A	N/A	N/A	N/A	N/A	N/A	N/A								N/A										N/A	N/A				
OtoTherm, shore 70		<input type="checkbox"/>	<input type="checkbox"/>	N/A	N/A	N/A	N/A	N/A	N/A	N/A								N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A			

Lacquering

Hard
 Soft
 Matt finish⁽⁴⁾
 Note: OtoTherm only in soft lacquering. Soft mould only in matt or soft.

Pull Out String

- Short - 5 mm
 Standard - 8 mm
 Long - 10 mm

Important Notes:

- ¹⁾ Not available for Mould and Micro Mould in hard acrylic material;
²⁾ Not available for Mould and Micro Mould;
³⁾ Black colour not available for Mould and Micro Mould in hard acrylic material;
⁴⁾ Not available for MicroShell

Default option. If customer doesn't choose any option then default will be chosen by production

Cosmetics & Usability

Comments